PORTUGUÊS

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Revisão: 13

PRODUCT

Technical Name: Box / Case for Dental Instruments

Commercial Name: Sterilization Autoclave Case

ATTENTION!

NON-STERILE PRODUCT.

The Autoclave Cases are supplied NON-STERILE and must be sterilized before use.

COMMERCIAL MODELS:

See Annex I: Commercial Models

1. INDICATION OF USE AND PURPOSE

The Sterilization Autoclavable Cases consist of boxes produced in polymer, indicated for conditioning medical products (Implacil instruments and implants) for sterilization. The use of cases optimizes the storage and organization of the instruments during surgery, contributing to its success and safety.

2. HOW TO USE

The Sterilization Autoclavable Cases must be sanitized before the sterilization process and its consequent inclusion in a sterilized surgical field.

Follow the instructions below:

CLEANING

Cleaning must be conducted in such a way that all parts of the item are exposed to the cleaning process.

To clean Implacil cases use only the following manual cleaning process.

Automated cleaning has not been validated. Do not use automated cleaning for Implacil cases.

1. Initial rinse with cold water. The initial cold rinse is carried out to remove any gross contamination present in the item. Cold water is used, as temperatures above 60°C / 140 ° F will clot proteins, making it difficult to remove them from contaminated items.

2. Prepare the enzymatic detergent according to the recommendation *.

3. Dip the products in the prepared detergent solution and let it soak for at least 5 minutes. While soaking, brush the items with a soft bristle brush.

4. Rinse: after removing the cases from the detergent solution, rinse thoroughly for at least 3 times (for at least 1 minute) under running water.

5. Visually inspect each product for visible dirt. If any visible dirt remains, repeat the cleaning process (steps 3 and 4) until there is no visible dirt left.

* Recommendation for the cleaning step

Always use distilled, deionized or demineralized water for this procedure.

Never use saline solutions, especially sodium hypochlorite and saline, disinfectants, hydrogen peroxide or alcohol for cleaning or rinsing the cutters.

Never use steel wool or sponges and abrasive products, so that the cutters are not damaged.

Cleaning: Enzyme Detergent Composition: Protease, Lipase and Amylase enzymes

Dilution: 1mL for each liter of water;

Action time: 5 min (Manual);

Water temperature: 40 ° C to 55 ° C (Manual);

Note: The cleaning process must be conducted in such a way that all parts of the cases are exposed. Contamination remaining on an item after cleaning reduces the effectiveness of any subsequent decontamination or sterilization process.

Clean the cases as soon as possible after use. Do not allow blood and waste to dry on cases. If cleaning needs to be delayed, place the case in a container covered with enzyme solution or appropriate detergent to delay drying.

The cases should be replaced as soon as the professional realizes their reduced effectiveness or degradation, according to the manufacturer of the raw material, the cases withstand more than 1,000 cycles of steam sterilization at 134 °C without any significant loss of properties.

STERILIZATION

This product must be properly sanitized and sterilized before use. Please use only steam sterilization according to the following parameters for sterilization:

Sterilization Parameters	
Temperature	132°C
Sterilization time	20 min
Drying time	20 - 60 min

3. COMMERCIAL PRESENTATION FORMS OF THE PRODUCT

The cases are packaged in a cardboard box and contain 1 unit of the product.

This instruction for use will be available in electronic format on the Implacil website (<https://implacilosstem.com.br/>). The printed instruction for use will be available to the customer who requests it and will be sent free of charge. The printed instruction for use can be requested through the Customer Service – 0800.779.0011, and will be sent at no additional cost. To locate it, check the “IFU” number and the version of the instruction that is indicated on the label of the product purchased. To read the instructions for use on the website, the Adobe Acrobat Reader program is required. User Alert: Always observe the correlation of the version of the instructions for use indicated on the label of the product purchased, as made available by Implacil.

4. PHYSICAL PRINCIPLE AND FUNDAMENTALS OF THE PRODUCT'S TECHNOLOGY, APPLIED TO ITS OPERATION AND ACTION:

The Sterilization Autoclavable Case consists of boxes produced in Polysulfone or Radel®, which offers excellent heat resistance and high mechanical resistance, indicated for conditioning medical products (Implacil instruments and implants) for sterilization. The use of cases optimizes the storage and organization of the instruments during surgery, contributing to its success and safety.

5. COMPOSITION

The following are used as raw material for the cases:

Radel® or Polysulfone.

Both raw materials have the following characteristics: translucent and ductile, high resistance and rigidity, high heat resistance (HDT), biocompatible and compatible with the sterilization method used.

6. EXPIRY DATE

5 years from the date of manufacture.

7. SPECIAL CONDITIONS FOR TRANSPORTATION, STORAGE AND HANDLING OF THE PRODUCT:

Storage: Keep it in the original packaging protected from rain and direct sunlight.

Transport: Protect the packaging from direct rainfall when transporting.

Packing: Keep in a place protected from rain and direct sunlight, at room temperature.

8. WARNINGS AND/OR PRECAUTIONS TO BE ADOPTED

The customer is fully responsible for the validation of the sterilization process. We recommend only the steam autoclave process for sterilization of the Sterilization Autoclavable Cases.

Standard applicable to sterilization in steam autoclave: *BS EN ISO 17665-1:2006 - Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices.*

We do not recommend sterilization by ethylene oxide, gamma radiation or dry heat. We recommend the immediate use after the sterilization process.

9. UNDESIRABLE SECONDARY OR SIDE EFFECTS AND CONTRAINDICATIONS:

The incorrect use or manipulation may cause damage and possible malfunction of the product.

10. PRECAUTIONS IN CASE OF DISPOSAL

To avoid environmental contamination or improper use of the Sterilization Autoclavable Cases, when these are no longer usable, they should be disposed in the appropriate

location, i.e., in a place for potentially contaminated materials.

11. WARRANTY TERMS

IMPLACIL OSSTEM - MATERIAL ODONTOLÓGICO S.A. guarantees the owner of this product warranty against any material or manufacturing defect that may appear within 365 days, counting from the original date of purchase. The product may have indefinite validity if transported, stored and packed suitably according to these Instructions for Use.

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validation and routine control of a sterilization process for medical devices.

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




22469 – ESTOJO CIRÚRGICO (UPGRADE)

34593 - ESTOJO STRETTO



Technical drawing of a dental implant assembly, showing the implant body, abutment, and crown. The drawing includes dimensions for the implant body (8mm, 10mm, 11mm, 13mm, 15mm) and the abutment (8mm, 10mm, 11mm, 13mm, 15mm). The crown is labeled "Fresa Escarada" and "Fresa Larga helicoidal". The assembly is labeled "Implacil OSSTEM" and "veloce".

The diagram illustrates the Primo Implantsystem components and their dimensions. The components are arranged in a row, each with a corresponding dimension list:

- Fresa Larga Helicoidal**: Dimensions: 8, 7, 10, 11, 12, 13, 14, 15.
- Tonel Check Core Cutlure**: Dimensions: 7, 10, 11, 12, 13, 14, 15.
- Fresa Conica 3.0**: Dimensions: 7, 10, 11, 12, 13, 14, 15.
- Fresa Conica**: Dimensions: 7, 10, 11, 12, 13, 14, 15.
- Excavador Ovoides/Conicofuglio**: Dimensions: 7, 10, 11, 12, 13, 14, 15.
- Tonel Check Conico**: Dimensions: 7, 10, 11, 12, 13, 14, 15.

The dimensions are listed in a column next to each component, indicating the length of the tool in millimeters.

The diagram illustrates the Implanon OSSTEM system, which consists of two main components: a top view and a bottom view.

Top View: This view shows a dark, rectangular device with rounded corners. It features a grid of 16 small white dots arranged in 4 rows and 4 columns. The Implanon OSSTEM logo is prominently displayed in the upper left corner. On the right side, the text "LEXANTAMBITA CEEBO" and "SOLUTRA" is visible. Two white, rectangular tabs are located on the left and right edges of the device.

Bottom View: This view shows the underside of the device, which is white with a grid of 16 circular holes arranged in 4 rows and 4 columns. Two vertical strips of dark, segmented material are positioned on the left and right sides of the device, likely serving as hinges or supports. The overall design is compact and rectangular.

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