**PORTUGUÊS**

Implante Universal III

**INGLÊS**

Implant Universal III

**ESPAÑHOL**

Implante Universal III

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Revisão: 08

## PRODUCT:

Technical name: Dental Implants (Osseointegrated)

Commercial name: Implant Universal III

## ATTENTION!!!

**SINGLE USE PRODUCT. REPROCESSING PROHIBITED.**

**DO NOT USE IF PACKAGING IS DAMAGED OR BROKEN.**

**STERILE PRODUCT. DO NOT RESTERILIZE. STERILIZED BY GAMMA RADIATION.**

**MANUFACTURER RECOMMENDS SINGLE USE.**

## Commercial models:

UN III Due Cone Maestro Implant  
UN III HI Maestro Tapered Implant  
UN III HE Maestro Tapered Implant  
UN III HI Tapered Implant  
UN III HE Tapered Implant  
UN III Morse Taper Implant  
UN III Morse Taper AR Due Cone Implant  
UN III Morse Taper CAR Implant  
UN III Switch Cylindrical HI Implant  
UN III Switch Cylindrical HE Implant  
UN III Switch Tapered HI Implant  
UN III Switch Tapered HE Implant  
UN III New Switch Cylindrical HI Implant  
UN III New Switch Cylindrical HE Implant  
UN III New Switch Tapered HI Implant  
UN III New Switch Tapered HE Implant  
UN III HI Implant  
UN III HE Implant  
UN III SC Morse Taper Implant  
UN III SC Morse Taper AR Implant  
UN III SC Morse Taper CAR Implant

## 1. INDICATION OF USE AND PURPOSE

Implacil implants are made from titanium in accordance with ASTM F-67 standard and are available in a variety of sizes.

Indication: Patients who require total or partial oral rehabilitation.

Purpose: To perform the function of the root of a missing tooth.

## 2. HOW TO USE THE PRODUCT

### 2.1 UN III Due Cone Maestro Implant

Implant

Cover

The UN III DUE CONE MAESTRO IMPLANT consists of a conical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. Further improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residues, causing the implant placement torque, inducing improvements in integration time. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Taper connection of this implant favors its use for unitary elements. It has an internal 2.00 mm hexagon used both for insertion and as an index for the protective elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of the Maestro Implant allows this implant to be installed in any type of bone with low insertion torque. Its healing chambers provide faster bone neoformation, favoring early loading of this implant, provided that the necessary prerequisites for this technique for Type I, II and III bones are followed, or delayed loading for any type of bone. Immediate loading is not indicated as standard, but its application is at the dentist's discretion, considering the clinical viability of the case. It is produced from turning titanium bars in accordance with the ASTM F-67 standard.

## Features:

Early or late rehabilitation;

Single/multiple;

Installation in any bone density: type I, II, III and IV;

Low bone tissue stitching during implant insertion;

Increases the osteotomy diameter;

Improves the quality of newly formed bone tissue;

Double Morse isolation that allows greater stability and greater bacterial isolation of the prosthetic component;

Includes 1 mm cover;

• Installation: Hex Key #7 – 1.17 mm;

• Implant Cover Placement Torque: 10Ncm

## Implant Application:

• Implant installation 2 mm below the bone;

• Gingival profile above 1.5/2 mm required;

## Rotation:

• Drilling: 600 rpm;

• Installation: 20 rpm;

## Suggested installation torque:

• 25 Ncm for Early\* / Late loading.

\*Early loading should be avoided in patients with:

- 1) Systemic diseases, such as diabetes, osteopenia, osteoporosis or those that cause changes in bone metabolism;
- 2) Occlusal dysfunction;
- 3) Implants installed in Type IV bone and/or in grafted areas.

## 2.2 UN III HI Maestro Tapered Implant

The UN III HI MAESTRO TAPERED IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The body of the implant is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. Still aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. In the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed uniformly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The internal hexagonal connection makes this system one of the simplest to handle. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of the Maestro Implant allows this implant to be installed in any type of bone with low insertion torque. Its healing chambers provide faster bone neoformation, favoring early loading of this implant, provided that the necessary prerequisites for this technique for Type I, II and III bones are followed, or delayed loading for any type of bone. Immediate loading is not indicated as standard, but its application is at the dentist's discretion, considering the clinical viability of the case. It is produced by turning titanium bars according to ASTM F-67 standard.

### FEATURES

- Increases the diameter of the osteotomy.
- Improves the quality of the newly formed bone tissue;
- Indication for single cases and safety for rehabilitation of multiple implants;

Supplied with 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
  - Implant Cover Placement Torque: 10 Ncm
- Implant application:
- Implant installation at bone level;

Rotation:

- Drilling: 600 rpm;
- Installation: 20 rpm;

Suggested installation torque:

- 25 Ncm for Early\* / Late loading.

\*Early loading should be avoided in patients with:

- 1) Systemic diseases, such as diabetes, osteopenia, osteoporosis or those that cause changes in bone metabolism;
- 2) Occlusal dysfunction;
- 3) Implants installed in Type IV bone and/or in grafted areas.

## 2.3 UN III HE Maestro Tapered Implant

The UN III HE MAESTRO TAPERED IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness titanium finish to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. Also in the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of the Maestro Implant allows this implant to be installed in any type of bone with low insertion torque. Its healing chambers provide faster bone neoformation, favoring early loading of

this implant, provided that the necessary prerequisites for this technique for Type I, II and III bones are followed, or delayed loading for any type of bone. Immediate loading is not indicated as standard, but its application is at the dentist's discretion, considering the clinical viability of the case. It is produced from the turning of titanium bars according to the ASTM F-67 standard.

#### FEATURES

Increases the diameter of the osteotomy.  
Improves the quality of the newly formed bone tissue.  
Indication for single cases and safety for rehabilitation of multiple implants;

Includes 1 mm cover;  
• Installation: Hex Key No. 7 – 1.17 mm;  
• Implant Cover Placement Torque: 10Ncm

Implant application:  
• Installation of the implant at bone level;

Rotation:  
• Drilling: 600 rpm; • Installation: 20 rpm;

Suggested installation torque:  
• 25 Ncm for Early\* / Late loading.

\*Early loading should be avoided in patients with:

- 1) Systemic diseases, such as diabetes, osteopenia, osteoporosis or those that cause changes in bone metabolism;
- 2) Occlusal dysfunction;
- 3) Implants installed in Type IV bone and/or in grafted areas.

#### 2.4 UN III HI Tapered Implant

The UN III HI TAPERED IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. In the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a neck that is slightly

divergent. The internal hexagonal connection makes this system one of the simplest to handle. In the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

Includes 1 mm cover;  
• Installation: Hex Key #7 – 1.17 mm;  
• Implant Cover Placement Torque: 10Ncm

Rotation:  
• Drilling: 800-1,200 rpm;  
• Installation: 20 rpm;

Rotation for 5/6mm long implants:

• Drilling: 200-300 rpm;  
• Installation: 20 rpm;

• Suggested installation torque of up to 60 Ncm

#### 2.5 UN III HE Tapered Implant

The UN III HE TAPERED IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. Also in the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread, which allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, and is indicated mainly for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard. The UN III Conical HE IMPLANT combines the efficiency of the external hexagonal prosthetic fitting,

recognized worldwide, with the practicality provided by the conical design.

Includes 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;

Rotation for 5/6 mm long implants:

- Drilling: 200-300 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.6 UN III Morse Taper Implant

The UN III MORSE TAPER IMPLANT consists of a conical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. This surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Taper connection of this implant favors its use for unitary elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being indicated mainly for use in low-density bone. It is produced from the turning of titanium bars according to ASTM F-67 standard.

Comes with a 1 mm cover;

- Installation: Hex Key #7 – 1.17 mm;
- Implant Cover Placement Torque: 10Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

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The UN III MORSE TAPER AR DUE CONE IMPLANT consists of a conical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. This surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Taper connection of this implant favors its use for unitary elements. It has an internal 2.00 mm hexagon used both for insertion and as an index for the prosthetic elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the fixation of the prosthetic component by a screw. The design of this implant favors load distribution, and is indicated mainly for use in low-density bone. It is produced by turning titanium bars according to ASTM F-67 standard.

Includes 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.8 UN III Morse Taper CAR Implant

The UN III MORSE TAPER CAR IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. This

surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Cone connection of this implant favors its use for unitary elements. It has an internal 2.00 mm hexagon used both for insertion and as an index for the prosthetic elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution and is indicated mainly for use in low-density bone. It is produced by turning titanium bars according to ASTM F-67 standard.

Includes 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.9 UN III Switch Cylindrical HI Implant

The UN III SWITCH CYLINDRICAL HI IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant is fully treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface. The implant has three millings at the apex, distributed evenly (120°) to facilitate its insertion. Its seating platform has the following diameters: 5.0 mm, which corresponds to the 4.75 mm implant. It is produced with a slightly divergent neck and a self-threading conical apex. All diameters presented in this model have a 2.5 mm internal hexagon and M2.0 internal thread. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. It is produced from the turning of titanium bars according to ASTM F-67 standard.

Includes a 1 mm cover;

- Installation: Hex Key #7 – 1.17 mm;
- Implant Cover Placement Torque: 10Ncm

Rotation:

- Drilling: 800-1,200 rpm;

- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.10 UN III Switch Cylindrical HE Implant

The UN III SWITCH CYLINDRICAL HE IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant is fully treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface. The implant has three millings at the apex, distributed evenly (120°) to facilitate its insertion. Its seating platform has the following diameters: 5.00 mm, which corresponds to the 4.75 mm implant. It is produced with a slightly divergent neck and a self-threading conical apex. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. All diameters presented in this model have a 2.7 mm external hexagon and M2.0 internal thread. On the longitudinal axis of the upper part of the implant there is a central thread that allows the fixation of the prosthetic component by a screw. It is produced by turning titanium bars according to ASTM F-67 standard. The UN III Switch Cylindrical HE Implant has the prosthetic fitting that has been the leader in the choice of professionals worldwide for decades. Its neck, or counter-sink, adapts well to the cortical bone, efficiently distributing the occlusal load.

Includes 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.11 UN III Switch Tapered HI Implant

The UN III SWITCH TAPERED HI IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant is fully treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface. In the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The internal hexagonal connection makes this system one of the simplest to handle. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. All diameters presented in this model have a 2.5mm internal hexagon and M2.0 internal thread. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

Includes 1 mm cover;

- Installation: Hex Key #7 – 1.17 mm;
- Implant Cover Placement Torque: 10Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.12 UN III Switch Tapered HE Implant

The UN III SWITCH TAPERED HE IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant is fully treated (sandblasted), a process that increases its

contact surface and facilitates anchoring by increasing the osseointegration surface. In the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. All diameters presented in this model have a 2.7mm external hexagon and M2.0 internal thread. The design of this implant favors load distribution, and is indicated mainly for use in low-density bone or even for immediate loading. It is produced from the turning of titanium bars according to ASTM F-67 standard. The UN III SWITCH TAPERED HE IMPLANT combines the efficiency of the external hexagonal prosthetic fitting, recognized worldwide, with the practicality provided by the conical design.

Includes 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.13 UN III New Switch Cylindrical HI Implant

The UN III NEW SWITCH CYLINDRICAL HI IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant is fully treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface. The implant has three millings at the apex, distributed evenly (120°) to facilitate its insertion. Its seating platform has the following diameters: 4.00 mm, which corresponds to the 3.75 mm implant; 5.0 mm, which corresponds to the 4.75 mm implant. It is produced with a slightly divergent neck and a self-threading conical apex. All diameters presented in this model have a 2.3 mm internal hexagon and M1.8 internal thread. On the longitudinal axis of the upper part of the implant there is a central thread that allows the fixation of

the prosthetic component by a screw. It is produced by turning titanium bars according to ASTM F-67 standard.

Includes 1 mm cover;

- Installation: Hex Key #7 – 1.17 mm;
- Implant Cover Placement Torque: 10Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

#### 2.14 UN III New Switch Cylindrical HE Implant

Implante

Cover

The UN III NEW SWITCH CYLINDRICAL HE IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant is fully treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface. The implant has three millings at the apex, distributed evenly (120°) to facilitate its insertion. Its seating platform has the following diameters: 4.00mm, which corresponds to the 3.75mm and 4.00mm implant; 5.00mm, which corresponds to the 4.75mm implant. It is produced with a slightly divergent neck and a self-threading conical apex. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. All diameters presented in this model have a 2.4mm external hexagon and M1.8 internal thread. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. It is produced by turning titanium bars according to ASTM F-67 standard. The UN III NEW SWITCH CYLINDRICAL HE IMPLANT has the prosthetic fitting that has been the leading choice of professionals worldwide for decades. Its neck, or counter-sink, adapts well to the cortical bone, efficiently distributing the occlusal load.

Includes 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

#### 2.15 UN III New Switch Tapered HI Implant

Implante

Tapa/Cover

The UN III NEW SWITCH TAPERED HI IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant is fully treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface. In the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The internal hexagonal connection makes this system one of the simplest to handle. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. All diameters presented in this model have a 2.3mm internal hexagon and M1.8 internal thread. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

Includes 1 mm cover;

- Installation: Hex Key #7 – 1.17 mm;
- Implant Cover Placement Torque: 10Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

#### 2.16 UN III New Switch Tapered HE Implant

Implante

Cover

The UN III NEW SWITCH TAPERED HE IMPLANT consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant is fully treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the

osseointegration surface. In the upper part, this implant has micro-spirals that improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. It is produced with a slightly divergent neck. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. All diameters presented in this model have a 2.4mm external hexagon and M1.8 internal thread. The design of this implant favors load distribution, and is indicated mainly for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard. The UN III NEW SWITCH TAPERED HE IMPLANT combines the efficiency of the external hexagonal prosthetic fitting, recognized worldwide, with the practicality provided by the conical design.

Includes 1 mm cover;

- Installation: Hex Key No. 7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;

- Suggested installation torque of up to 60 Ncm

### 2.17 UN III HI Implant

The UN III HI IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. The implant has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. Its seating platform has the following diameters: 3.50 mm, which corresponds to the 3.3 mm implant; 4.00 mm, which corresponds to the 3.75 mm and 4.3 mm implant; 5.0 mm, which corresponds to the 4.75 mm implant. It is produced with a slightly divergent neck and a self-threading conical apex. The internal hexagonal connection makes this system one of the simplest to handle. On the longitudinal axis of the upper part of the implant

there is a central thread that allows the prosthetic component to be fixed by a screw. It is produced by turning titanium bars according to ASTM F-67 standard.

Includes 1 mm cover;

- Installation: No. 7 Hex Key – 1.17 mm;
- Implant Cover Placement Torque: 10Ncm

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;

- Suggested installation torque of up to 60 Ncm

### 2.18 UN III HE Implant

The UN III HE IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. The implant has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. Its seating platform has the following diameters: 3.50 mm, which corresponds to the 3.3 mm implant; 4.00 mm, which corresponds to the 3.75 mm and 4.30 mm implant; 5.00 mm, which corresponds to the 4.75 mm implant. It is produced with a slightly divergent neck and a self-threading conical apex. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. It is produced by turning titanium bars according to ASTM F-67 standard. The UN III Cylindrical HE IMPLANT has the prosthetic fitting that has been the leading choice of professionals worldwide for decades. Its neck, or counter-sink, adapts well to the cortical bone, efficiently distributing the occlusal load.

Includes 1 mm cover;

- Installation: Hex Key #7 – 1.17 mm;
- Implant Cover Placement Torque: 10 Ncm

Rotation:

- Drilling: 800-1,200 rpm;

- Installation: 20 rpm;
- Suggested installation torque of up to 60 Ncm

### 2.19 UN III SC Morse Taper Implant

The UN III SC MORSE TAPER IMPLANT consists of a conical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. This surface treatment is carried out on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Taper connection of this implant favors its use for unitary elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the fixation of the prosthetic component by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone. It is produced from the turning of titanium bars according to ASTM F-67 standard. This implant model does not come with an implant cover.

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Installation: Hex Key #7 – 1.17 mm;
- Suggested installation torque of up to 60 Ncm

### 2.20 UN III SC Morse Taper AR Implant

The UN III SC MORSE TAPER AR IMPLANT consists of a conical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. This surface treatment is performed on the entire implant, since

it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Taper connection of this implant favors its use for unitary elements. It has an internal 2.00 mm hexagon used both for insertion and as an index for the prosthetic elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the fixation of the prosthetic component by a screw. The design of this implant favors load distribution, and is mainly indicated for use in low-density bone. It is produced by turning titanium bars according to ASTM F-67 standard.

This implant model does not come with an Implant Cover.

Rotation:

- Drilling: 800-1,200 rpm;
- Installation: 20 rpm;
- Installation: Hex Key #7 – 1.17 mm;
- Suggested installation torque of up to 60 Ncm

### 2.21 UN III SC Morse Taper CAR Implant

The UN III SC MORSE TAPER CAR IMPLANT consists of a cylindrical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. This surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Cone connection of this implant favors its use for unitary elements. It has an internal 2.00 mm hexagon used both for insertion and as an index for the prosthetic elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the fixation of the prosthetic component by a screw. The design of this implant favors load distribution, and is mainly indicated for use in low-density bone. It is produced by turning titanium bars according to ASTM F-67 standard.

This implant model does not come with an Implant Cover.

Rotation:

- Drilling: 800-1,200 rpm;

- Installation: 20 rpm;
- Installation: Hex Key #7 – 1.17 mm;
- Suggested installation torque of up to 60 Ncm

Note: See Annex I for dimension references; Annex II for the compatibility table and Annex III for the milling sequences.

### 3. PRODUCT PRESENTATION FORM

The product is presented in individual packaging.

Package contents:

- 01 Implant Universal III (UN III Due Cone Maestro Implant, UN III HI Maestro Tapered Implant, UN III HE Maestro Tapered Implant, UN III HI Tapered Implant, UN III HE Tapered Implant, UN III Morse Taper Implant, UN III Morse Taper AR Due Cone Implant, UN III Morse Taper CAR Implant, UN III Switch Cylindrical HI Implant, UN III Switch Cylindrical HE Implant, UN III Switch Tapered HI Implant, UN III Switch Tapered HE Implant, UN III New Switch Cylindrical HI Implant, UN III New Switch Cylindrical HE Implant, UN III New Switch Tapered HI Implant, UN III New Switch Tapered HE Implant, UN III HI Implant, UN III HE Implant, UN III SC Morse Taper Implant, UN III SC Morse Taper AR Implant, UN III SC Morse Taper CAR Implant).

- 01 Implant Cover (OPTIONAL)
- 03 Traceability Labels
- 01 Assembler (OPTIONAL)
- 01 Assembler Screw (OPTIONAL)
- 01 Certificate of Origin

These instructions for use will be available on the Implacil website (<https://implacilosstem.com.br/>) in accordance with current legislation on the provision of instructions for use in non-printed format for health products – RDC 751/2022. Printed instructions for use will be available to customers who request them and will be sent free of charge.

To locate it, check the “IFU” number and the version of the instruction that is indicated on the Label of the purchased product.

To read the instructions for use on the website, you need Adobe Acrobat Reader.

User Alert: Always observe the correlation of the version of the instructions for use indicated on the label of the purchased product, as made available by Implacil.

The certificate of origin is sent for traceability control by the patient – the dentist must place one of the traceability labels on the certificate and give it to the patient.

### 4. PHYSICAL PRINCIPLE AND FUNDAMENTALS OF THE PRODUCT'S TECHNOLOGY, APPLIED TO ITS OPERATION AND ACTION

The operating principle of a dental implant is the biocompatibility of the material from which it is made, which allows its bone integration (maxillary and mandibular), performing the function of the root of the lost

tooth. The thread is designed to conveniently distribute the masticatory load. The porosity of the surface exerts bone retention, which is distributed homogeneously in the surface area and also acts in the micro distribution of the masticatory load. The apical recesses act after osseointegration as mechanical retention.

The implant is inserted into the bone through an intraoral surgical process, fixed by locking in bone type I, II, III or IV. Implacil implants were designed for single, multiple, total restorations and also for the Immediate Loading system. Check the application indication for each model in item 2.

### 5. OPTIONALS, CONSUMABLES AND SUPPORT MATERIALS

**Optional:** Not applicable

**Consumables:** not applicable

**Support materials:** Instructions for use – available electronically at: <https://implacilosstem.com.br/>

**Manufacturing items for exclusive use in the company's products:**

All parts described in these Instructions for use are for exclusive use in the IMPLANTE UNIVERSAL III product.

Attention: The use of any part, accessory or material not specified or provided for in these Instructions for use is the sole responsibility of the user.

### 6. INSTRUCTIONS FOR USE OF THE PRODUCT

- Open the outer packaging. Check that the blister is perfectly sealed. Reject the product if the packaging is damaged.
- Start the sequence of the drills (see Annexe III) to create the cavity where the implant will be inserted. Open the blister only at the time of insertion of the implant.
- Remove the capsule from the blister and open it carefully. Do not throw away the cap, as the implant cap is inside it.
- Remove the implant from inside the capsule using the key.
- Insert the implant into the newly drilled cavity. Maximum insertion speed of 20 RPM. If done manually, the maximum torque to be used is 60 N.mm.
- Once inserted, remove the mounter (when applicable) using a 1.17mm hex key.
- Remove the implant cap from the capsule lid.
- Seal with the implant cover.
- Installation complete.

### 7. COMPOSITION

As raw material for implants, it is produced from the turning of titanium bars – ASTM F67.

### 8. VALIDITY PERIOD

05 years from the date of sterilization.

### 9. SPECIAL CONDITIONS FOR TRANSPORTATION, STORAGE AND PACKAGING OF THE PRODUCT:

Store at room temperature, in a cool place, away from light and heat. Avoid exposure to heat and excessive humidity.

Do not place excessive load on the packaging to avoid compromising the seals that guarantee the sterility of the product.

The Implant Universal III Implacil must be transported in a clean, dry place at room temperature. Transportation must be carried out appropriately to avoid falling and damaging the original packaging and mechanical shocks that may damage the product.

**ATTENTION !!!**

**Single-use product. Reprocessing prohibited.**

**Do not use if packaging is broken or damaged.**

**Sterile Product. Do not resterilize. Sterilized by gamma radiation.**

**Manufacturer recommends single use.**

**10. WARNINGS AND/OR PRECAUTIONS TO BE TAKEN**

For the safe and effective use of Implacil Universal III Implants, we strongly recommend that professionals undergo specific training (qualified professional). The surgical techniques required to place these implants are highly specialized and complex. The use of an improper technique may result in implant failure, loss of supporting bone or other adverse effects.

**Performance on essential safety and efficacy requirements of the medical product, and possible undesirable side effects:**

**10.1 PATIENT SCREENING AND EVALUATION**

- Rigorous screening of implant candidates should consider factors such as age and parafunctional habits and physiological and anatomical conditions.
- The patient's general health status is essential, with special attention to allergies to medications, local or systemic conditions that may interfere with the healing or osseointegration processes.

**10.2 SELECTION OF IMPLANTS**

The choice of model, diameter, length, position and quantity of implants must be made individually, considering:

- Anatomy and region to be rehabilitated;
- Bone quality and quantity (adequate bone thickness for implant installation);
- Available space.

**10.3 IMPORTANCE OF IMAGING EXAMS**

Radiographic exams are essential to evaluate the patient's bone and dental structure.

Types of recommended radiographs:

- Panoramic radiograph
- Periapical radiograph

- Computed tomography

Purposes:

- Identify injuries, infections, diseases or fractures;
- Determine bone quality and quantity;
- Locate important anatomical structures;
- Evaluate proximity of adjacent teeth;
- Define exact incision points;
- Detect corrosion, gingival and bone damage.

Note: X-rays are safe after implant installation, as long as they are performed as recommended by your dentist.

**10.4 RECOMMENDATIONS FOR USE**

Sterilization and Reuse:

- Single-use product; should not be reprocessed, reused or resterilized. Reuse may cause biological contamination.
- Opened or damaged packaging should not be used.
- Make sure the product is within its expiration date.

Surgical Technique:

- The critical threshold for bone necrosis is 47°C for one minute; above 50°C for more than one minute, irreversible loss of bone regeneration capacity occurs.
- Always use the Implacil product sequence for implant installation and activation.
- Before each procedure, check the fit of the parts and the condition of the instruments. Replace worn ones.
- Use cutting burs, with progressive diameter, with appropriate speed and abundant irrigation to avoid overheating.
- For short ST implants (5 mm/6 mm), assess the implant-crown relationship; the system does not allow angled components.
- Check the recommended installation torque. In case of high torque, the driver may get stuck inside the implant; a slight counter-torque may solve this.
- Make sure the motor is calibrated to apply the correct torque.
- In cases of bone graft, the implant must be in contact with at least 1/3 of the patient's bone.

Post-Surgical Monitoring:

- The patient must be monitored after implant installation.
- Lack of primary stability may compromise osseointegration.
- Failures in surgical and/or prosthetic planning may lead to loss of the implant, fracture of components or loosening of prosthetic screws.

## 10.5 PRODUCT SAFETY AND EFFICACY

The Universal Implant III is completely safe, as long as safety rules are followed in all procedures performed. The product was developed with non-toxic medical-grade material, manufactured in titanium according to ASTM F-67 standard. When maintained (stored) and preserved as mentioned in item 9 of these Instructions for Use, the product will not lose or change its physical and dimensional characteristics. The product is sterile and suitable for use for up to 5 years from the date of sterilization.

### MAGNETIC RESONANCE IMAGING (MRI)

Implacil has not evaluated the safety of its dental implants in the MRI environment. The safety of magnetically induced displacement force and the extent to which radiofrequency (RF) induced heating is capable of affecting patient safety have not been evaluated. The safety of Implacil dental implants in the MRI environment is unknown. Performing an MRI on a patient with this device may result in harm to the patient.

Implacil recommends that medical professionals and the manufacturer of the MRI equipment be consulted for compatibility with implants prior to use. We inform you that Implacil implants are machined in ASTM F67 Gr.4 Titanium.

## 11. SIDE EFFECTS OR UNDESIREBLE EFFECTS AND CONTRAINDICATIONS

The implant installation surgery should only be performed by professionals who have adequate technical training, as it is a highly specialized and complex surgical technique.

This product is contraindicated for patients who show signs of allergy or hypersensitivity to the chemical components of the material: titanium.

Excessive bone loss or implant fracture may occur when it is loaded beyond its functional capacity.

Physiological and anatomical conditions may impair the performance of dental implants. This must be taken into consideration when dental implants are placed in patients with the following conditions:

Conditions that may Impair Implant Performance:

- Acute inflammatory or infectious processes in living tissue;
- Low-quality bone or bone insufficiency;
- Remains of roots at the implant insertion site;

- Poor-quality oral hygiene;
- Metabolic or systemic diseases that affect bone healing (diabetes, hypertension, coagulation disorders);
- Chemotherapy or radiotherapy treatments;
- Occlusion or joint problems;
- Uncontrolled periodontitis;
- Bruxism;
- Xerostomia;
- Excessive use of alcohol or tobacco;
- Pregnancy.

Possible Adverse Events:

- Implant integration failure or loss;
- Dehiscence requiring bone graft;
- Perforation of anatomical structures (maxillary sinus, inferior alveolar canal);
- Infection (abscess, fistula, inflammation);
- Persistent pain, paresthesia;
- Excessive bone loss requiring intervention;
- Implant fracture.

Failure of osseointegration may be caused by factors such as inadequate osteotomy, infections, adverse systemic conditions, compromised blood supply, inadequate instruments, poor oral hygiene, occlusal trauma or lack of specific training.

The success of dental implants depends on adequate planning, correct technique and careful post-operative follow-up.

## 12. INSTALLATION OR CONNECTION TO OTHER PRODUCTS

The UNIVERSAL III IMPLANT must be connected to a specific driver to be implanted. After being integrated, it receives the connection of prosthetic components so that the patient's oral restoration can be performed.

Attention: The use of any part, accessory or material not specified or provided for in these Instructions for Use is the sole responsibility of the user.

Implacil implants are only compatible with the components of our system.

Note: Consult the compatibility table in Appendix II.

## 13. PRECAUTIONS IN CASE OF DISPOSAL

To avoid environmental contamination or improper use of the IMPLANT UNIVERSAL III, when it is no longer in use, it must be discarded in an appropriate place,

that is, in a place for potentially contaminated materials.

**14. WARRANTY TERMS**

IMPLACIL OSSTEM - MATERIAL ODONTOLOGICO SA guarantees the owner of this product a guarantee against any material or manufacturing defect that may arise within 365 days, counting from the date of purchase by the first purchaser.

**15. TERM OF RESPONSIBILITY**

The company IMPLACIL OSSTEM - MATERIAL ODONTOLOGICO S.A., represented by its technical manager Danielle Aparecida Mariano Aviz and its legal representative Aluizio Leal do Canto Junior, assume technical and legal responsibility for the product and that all information provided regarding the product IMPLANTE UNIVERSAL III, contained in these Instructions for Use, is true.

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