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IFU-038

Revisão: 04

PORTUGUÊS

Implante Maestro

INGLÊS

Implant Maestro

ESPAÑHOL

Implante Maestro

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Revisão: 04

PRODUCT

Technical Name: Dental Implant (Osseointegrated)

Commercial Name: Implant Maestro

ATTENTION!

Single-use product. Reprocessing prohibited.

Do not use if packaging is broken or damaged.

Sterile Product. Do not resterilize. Sterilized by gamma radiation.

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Implant Maestro HI Switch

1. INDICATION OF USE AND PURPOSE

Implants are indicated for patients who require total or partial oral rehabilitation and are manufactured from titanium, in accordance with ASTM F-67 standard.

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The Implant Due Cone Maestro Superiore consists of a conical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. Also aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residues, reducing the implant placement torque, inducing an improvement in integration time. The Implant Due Cone Maestro Superiore has a differentiated surface treatment that involves a chemical surface treatment, a treatment that will improve the integration of the bone in the implant. This surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest.

The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. It is produced with a slightly divergent neck. The Morse Taper connection of this implant favors its use for single elements. It has an internal 2.00 mm hexagon used both for insertion and as an index for the prosthetic elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone. It is produced from the turning of titanium bars according to ASTM F-67 standard.

The Implant Due Cone Maestro Superiore is supplied with a 1 mm Healing Abutment, a component used in place of a cover screw to protect the prosthetic connection during the healing period and to promote soft tissue conditioning. This Healing Abutment is registered with ANVISA under a separate process, should the dental surgeon need to purchase it individually

2.2 Implant Maestro HE Superiore

The IMPLANT MAESTRO HE SUPERIORE consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. The "Superiore" line has a differentiated surface treatment that involves a chemical surface treatment, a treatment that will improve the integration of the bone in the implant. This surface treatment is applied to the entire implant, as it is recommended for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. It is produced with a slightly divergent neck. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be

fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

2.3 Implante Maestro HE Switch Superiore


The IMPLANTE MAESTRO HE SWITCH SUPERIORE consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that provide better integration; these channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. The “Superiore” line has a differentiated surface treatment that involves a chemical surface treatment, a treatment that will improve the integration of the bone in the implant. This surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

Platform Switch:

Developed with the objective of standardizing the prosthetic interface and facilitating the selection of abutments and other components by the professional, allowing the use of the same abutments on implants of different diameters.

This configuration occurs when the implant platform has a reduced diameter relative to the implant body, as shown in the example below:

2.4 Implante Maestro HI Superiore

The IMPLANTE MAESTRO HI SUPERIORE consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. The “Superiore” line has a differentiated surface treatment that involves a chemical surface treatment, a treatment that will improve the integration of the bone in the implant. This surface treatment is applied to the entire implant, as it is recommended for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. It is produced with a neck that is slightly divergent. The internal hexagonal connection makes this system one of the easiest to handle. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, and is mainly recommended for use in low-density bone or even for immediate loading. It is produced by turning titanium bars in accordance with ASTM F-67 standard.

2.5 Implante Maestro HI Switch Superiore


The IMPLANTE MAESTRO HI SWITCH SUPERIORE consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness finish of the titanium to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that provide better integration; these channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. The “Superiore” line has a differentiated surface treatment that involves a chemical surface treatment, a treatment that will improve the integration of the bone in the implant. This surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. The internal hexagonal connection makes this system one of the easiest to handle. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

Platform Switch:

Developed with the objective of standardizing the prosthetic interface and facilitating the selection of abutments and other components by the professional, allowing the use of the same abutments on implants of different diameters.

This configuration occurs when the implant platform has a reduced diameter relative to the implant body, as shown in the example below:

Exemplo HE

2.6 Implant Maestro HE Tapered

The Implant Maestro He Tapered consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness titanium finish to improve its adhesion to the bone tissue. Still aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residues, reducing the implant placement torque, inducing an improvement in integration time. Still in the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. It is produced with a slightly divergent neck. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

2.7 Implant Maestro HI Tapered

The Implant Maestro HI Tapered consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness titanium finish to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide better integration; such channels function as pockets of autogenous bone residues, reducing the implant placement torque, inducing an improvement in integration time. In the upper part, this implant has micro-spirals, which improve distribution in the bone cortex. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. It is produced with a slightly divergent neck. The internal hexagonal connection makes this system one of the easiest to handle. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

2.8 Implant Due Cone Maestro

The Implant Due Cone Maestro consists of a conical body with threads that aid its installation and anchoring in the bone bed along its entire length. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring, increasing the osseointegration surface. Also aiming at improving osseointegration, this model has radial channels that differentiate this implant profile from similar models that are already registered, as these channels provide

better integration; such channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. This surface treatment is performed on the entire implant, since it is indicated for installation 1 to 2 mm below the bone crest. The implant also has three millings at the apex, distributed evenly (120°) in order to facilitate its insertion. It is produced with a slightly divergent neck. The Cone Morse connection of this implant favors its use for single elements. It has an internal 2.00 mm hexagon used both for insertion and as an index for prosthetic elements. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, being mainly indicated for use in low-density bone. It is produced by turning titanium bars according to ASTM F-67 standard.

The Implant Due Cone Maestro is supplied with a 1 mm Healing Abutment, a component used in place of a cover screw to protect the prosthetic connection during the healing period and to promote soft tissue conditioning. This Healing Abutment is registered with ANVISA under a separate process, should the dental surgeon need to purchase it individually.

2.9 Implants Maestro HE Switch


The Implants Maestro HE Switch consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness titanium finish to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that provide better integration; these channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. The external hexagonal connection is a well-established prosthetic fitting, used worldwide for over 40 years. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution, and is mainly indicated for use in low-density bone or even for immediate loading. It is produced by

turning titanium bars in accordance with ASTM F-67 standard.

Platform Switch:

Developed with the objective of standardizing the prosthetic interface and facilitating the selection of abutments and other components by the professional, allowing the use of the same abutments on implants of different diameters.

This configuration occurs when the implant platform has a reduced diameter relative to the implant body, as shown in the example below:

2.10 Implante Maestro HI Switch


The Implante Maestro HI Switch consists of a conical body with threads along its entire length that aid its installation and anchoring in the bone bed. The implant body is mechanically treated (sandblasted), a process that increases its contact surface and facilitates anchoring by increasing the osseointegration surface, maintaining only in the neck region a low-roughness titanium finish to improve its adhesion to the bone tissue. Also aiming at improving osseointegration, this model has radial channels that provide better integration; these channels function as pockets of autogenous bone residue, reducing the implant placement torque, inducing an improvement in integration time. The implant also has three millings at the apex, distributed evenly (120°) to facilitate its insertion. The internal hexagonal connection makes this system one of the simplest to handle. On the longitudinal axis of the upper part of the implant there is a central thread that allows the prosthetic component to be fixed by a screw. The design of this implant favors load distribution and is mainly indicated for use in low-density bone or even for immediate loading. It is produced by turning titanium bars according to ASTM F-67 standard.

Platform Switch:

Developed with the objective of standardizing the prosthetic interface and facilitating the selection of abutments and other components by the professional, allowing the use of the same abutments on implants of different diameters.

This configuration occurs when the implant platform has a reduced diameter relative to the implant body, as shown in the example below:

Note: See Annex I for reference code and commercial model; Annex II for Compatibility; Annex III for Drill Sequence

3. FORMS OF COMMERCIAL PRESENTATION OF THE PRODUCT

The implants are packaged inside a capsule that is inserted into a blister sealed with Tyvek. The blister is placed inside a cardboard box.

The product is presented in individual packaging.

Package contents:

Model: Implant Due Cone Maestro Superiore, Implant Maestro HE Superiore, Implant Maestro HE Switch Superiore, Implant Maestro HI Superiore, Implant Maestro HI Switch Superiore, Implant Maestro HE Tapered, Implant Maestro HI Tapered, Implant Due Cone Maestro, Implant Maestro HE Switch, Implant Maestro HI Switch.

- 01 Implant
- 01 Implant Cover

All models come with 3 traceability labels, which must be attached: to the clinical record, to the document to be given to the patient, and to the tax documentation that generates the charge.

Each traceability label contains:

- a) commercial name or model;
- b) identification of the manufacturer;
- c) product or system component code;
- d) batch number and ANVISA registration number.

This instruction for use will be available on the Implacil website (<https://implacilosstem.com.br/>) in accordance with current legislation on the provision of instructions for use in non-printed format for health products – RDC 751/2022. Printed instructions for use can be requested through Customer Service – 0800.779.0011, and will be sent at no additional cost.

To locate it, check the "IFU" number and the version of the instruction that is indicated on the Label of the purchased product.

To read the instructions for use on the website, you need Adobe Acrobat Reader.

User Alert: Always observe the correlation of the version of the instructions for use indicated on the label of the purchased product, as made available by Implacil.

4. PHYSICAL PRINCIPLE AND FUNDAMENTALS OF THE PRODUCT'S TECHNOLOGY, APPLIED TO ITS OPERATION AND ACTION:

The operating principle of a dental implant is the biocompatibility of the material from which it is made, which allows its bone integration (maxillary and mandibular), performing the function of the root of the lost tooth. The thread is designed to conveniently distribute the masticatory load. The porosity of the surface exerts bone retention, which is distributed homogeneously in the surface area and also acts in the microdistribution of the masticatory load. The apical recesses act after osseointegration as mechanical retention.

The implant is inserted into the bone through an intraoral surgical process, fixed by locking in bone type I, II, III or IV. IMPLACIL implants were designed for single, multiple, total restorations and also for the Immediate Loading system.

5. OPTIONALS, CONSUMABLES AND SUPPORT MATERIALS

Optional: Not applicable

Consumables: Not applicable

Support materials: Instructions for Use

Manufacturing items for exclusive use in the company's products:

All parts described in these Instructions for Use are for exclusive use in the IMPLANTE product.

Attention: The use of any part, accessory or material not specified or provided for in these Instructions for Use is the sole responsibility of the user.

6. INSTRUCTIONS FOR USE OF THE PRODUCT

- a) Open the outer packaging. Check that the blister is perfectly sealed. Reject the product if the packaging is damaged.
- b) Start the sequence of the drills (see Annex III) to create the cavity where the implant will be inserted. Open the blister only at the time of insertion of the implant.
- c) Remove the capsule from the blister and open it carefully. Do not throw away the cap, as the implant cap is inside it.
- d) Remove the implant from inside the capsule using the key.
- e) Insert the implant into the newly drilled cavity. Maximum insertion speed of 20RPM. If done manually, the maximum torque to be used is 40N.cm
- f) Remove the implant cap from the capsule lid.

g) Seal with the implant cover.

h) Complete installation

7. COMPOSITION

The product was developed with non-toxic medical grade material, manufactured in titanium ASTM F-67.

8. VALIDITY PERIOD

05 years from the date of manufacture.

9. SPECIAL CONDITIONS FOR TRANSPORTATION, STORAGE AND PACKAGING OF THE PRODUCT:

Storage: Keep in a place protected from rain and direct sunlight, and in its original packaging.

Transport: When transporting, protect the packaging from direct rain.

Packaging: Keep in a place protected from rain and direct sunlight.

ATTENTION!!!

Single-use product. Reprocessing prohibited.

Do not use if packaging is broken or damaged.

Sterile Product. Do not resterilize. Sterilized by gamma radiation.

10. WARNINGS AND/OR PRECAUTIONS TO BE ADOPTED

For the safe and effective use of **Implacil implants**, we strongly suggest that specific training be undertaken, since the surgical techniques required to place these dental implants are highly specialized and complex procedures. Using an improper technique may lead to implant failure and/or loss of supporting bone or other adverse effects.

A rigorous screening of implant candidates should be performed, including age and the presence of parafunctional habits.

Visual inspection and periapical radiographs are essential for determining anatomical accidents, occlusion conditions, periodontal status and bone adequacy.

Cephalometric radiography, tomographies

- Do not sterilize or use a second time!

- Do not use an open or damaged package.

- The critical threshold (temperature/time relationship) for bone tissue necrosis is a temperature of 47°C applied for one minute. At a temperature of 50°C for more than one minute, the critical threshold for non-regeneration of peri-implant bone is reached.

- The use of cutting drills, with progressive diameter, at an appropriate speed (600 RPM/min), associated with abundant irrigation, allows excessive thermal elevation to be avoided.

11. PERFORMANCE OF ESSENTIAL REQUIREMENTS FOR SAFETY AND EFFICACY OF THE MEDICAL PRODUCT, AND ANY UNDESIRABLE SIDE EFFECTS

11.1 - INDICATION, PURPOSE OR USE FOR WHICH THE PRODUCT IS INTENDED

Indication: Patients who require total or partial oral rehabilitation.

Purpose: To perform the function of the root of the missing tooth

11.2 - UNDESIRABLE SIDE EFFECTS AND CONTRAINDICATIONS

Implant installation surgery should only be performed by professionals who have adequate technical training, as it is a highly specialized and complex surgical technique.

Excessive bone loss or implant fracture can occur when the implant is loaded beyond its functional capacity. Physiological and anatomical conditions can impair the performance of dental implants. This should be taken into consideration when placing dental implants in patients with the following conditions: poor bone quality; poor oral hygiene; medical conditions such as blood disorders or uncontrolled hormonal conditions; metabolic or systemic diseases that affect bone or wound healing, such as uncontrolled diabetes and coagulation disorders; chemotherapy or radiotherapy treatments; untreatable occlusion or joint problems; cases of insufficient bone height or width ; tobacco and alcohol abuse.

Potential adverse events associated with the use of dental implants may include: nonintegration; loss of integration; dehiscence requiring bone grafting; perforation of the maxillary sinus, inferior border, lingual surface or inferior alveolar canal; infection presented by: abscess, fistula, suppuration, inflammation or radiolucency; persistent pain, numbness, paresthesia; excessive bone loss requiring intervention; and implant fracture.

11.3 - PRODUCT SAFETY AND EFFICACY:

The Implant is completely safe, as long as safety rules are followed in all procedures performed. The product was developed with non-toxic medical-grade material, manufactured in titanium ASTM F-67. When maintained (stored) and preserved as mentioned in these Instructions for Use, the product will not lose or alter its physical and dimensional characteristics. The product is sterile and suitable for use for up to 5 years from the date of sterilization.

Magnetic Resonance Imaging (MRI): Implacil has not evaluated the safety of its dental implants in the MRI environment. The safety of magnetically induced displacement forces and the extent to which radiofrequency (RF) induced heating may affect patient safety have not been evaluated. The safety of Implacil dental implants in the MRI environment is unknown. Performing an MRI on a patient with this device may result in harm to the patient.

Implacil recommends that medical professionals and the manufacturer of the MRI equipment be consulted for compatibility with implants prior to use. Please be advised that Implacil implants are machined in Titanium Gr.4 ASTM F67.

12. INSTALLATION OR CONNECTION TO OTHER PRODUCTS

The IMPLANT, to be implanted, needs to be connected to a specific key, after being integrated it receives the connection of prosthetic components so that the patient's oral restoration can be carried out.

Attention: The use of any part, accessory or material not specified or provided for in these Instructions for Use is the sole responsibility of the user.

13. PRECAUTIONS IN CASE OF DISPOSAL

To avoid environmental contamination or misuse of the Implant, when it is no longer usable, it must be disposed of in hospital waste or in a location for potentially contaminated materials.

14. WARRANTY TERMS

IMPLACIL OSSTEM - MATERIAL ODONTOLOGICO SA guarantees the owner of this product against any material or manufacturing defect that may arise within 365 days, counting from the date of purchase by the first purchaser. It is estimated that dental implants have a useful life of approximately 15 years, therefore, as determined by Law No. 8078 of September 11, 1990, as long as there is a complaint within the ninety-day statute of limitations after the defect is evidenced, even if the defect only manifests itself after the end of the contractual warranty period, Implacil undertakes to replace the device with an identical or similar device.

15. TERM OF RESPONSABILITY

The company IMPLACIL OSSTEM - MATERIAL ODONTOLOGICO S.A., represented by its technical manager Danielle Aparecida Mariano Aviz and its legal representative Aluizio Leal do Canto Junior, assume technical and legal responsibility for the product and that all information provided regarding the IMPLANT MAESTRO product, contained in these Instructions for Use, is true.

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EN Note²: There is compatibility for the Superiore line, but restricted to the diameters previously determined for the respective models. The same applies to the following tables.

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MA= Moldeira Aberta / Open Tray / Bandeja abierta

MF= Moldeira Fechada / Closed Tray / Bandeja cerrada

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EN:

Implacil certifies that the connection system of the Implants that make up this family is safe and their combination with other products or equipment as instructed in the instructions for use does not alter the expected performance.

Implants from this family can only be used in interface with products manufactured by Implacil , whose interface codes are included in the product's Instruction for Use.

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